



# Sound Reducing Enclosure

## SONIFIER Acoustic Enclosure

Cat. # 20268

Although ultrasound is above the audible range of the human ear, mechanical noise occurs when liquids are treated ultrasonically. The noise may be disturbing in a room with low ambient noise level if exposed for prolonged periods of time.

Suitable for use with standard probes\* controlled via a Branson 250 or 450 Sonifier, a *Consonic-designed and manufactured* Sound Reducing Enclosure can minimize mechanical noise produced by ultrasonic processing.



CLEAR FRONT DOOR

- powder-coated metal cabinet
- lined with waterproof sound-abating material
- Easily cleaned
- Fully transparent door
- opening at the top enables the use of standard probes

*\* not for use with cup horns*

**Consonic**

WWW.CONSONIC.COM